Cap-Chur

Table of ContentsAlphabetized Listing

	•
A	
Arrow Adaptor	
B Blank Adaptor.22	
C	
Carrying Case - Pistol	
Carrying Case - Rifle	
Charges	
CO2 Gas Packs	
G	
Gas Packs	
Guns and Accessories	43-
M Marking Syringe Complete	
N	
Needles & Syringes	45-
Noseplug Barb	
Noseplug Heavy	
Noseplug Large Collar	
Noseplug Plain	
Noseplug Small Collar	
Noseplugs P	
Pistol	
Plunger Lube	
Positioner	
Powerloads	
Powerloads - Low	
Powerloads - Medium	
Powerloads - Very Low	
Powerloads & Charges	
Practice Syringe	
Projectors	
Projectors - Accessories	43-
R Rifle - Long Range	
Rifle - X-Long Range	
Rubber Plunger	
S	
Sideport	45-
Sideport Noseplug - Drop Out	
Sideport Small Collar	
Stock Extension	
Syringe Barrel Only	
Syringes	
T Tail Piece	

Guns an	d Accessories	5	
Projecto	rs		
Pistol	Short range projector: a pi	istol type projector which is powered by one CO2 gas pack d range projector: a pistol type projector that is breech load	
	KVS 039-305 039-307	Size Short Range Mid Range	Mfg. Cap-Chur Cap-Chur
Rifle - Long	Long range projector: a rif	le type projector which is powered by two CO2 gas packs. ⁻ yringes a distance of 25-30 yards with improved accuracy o Size Long Range	
Rifle - X-Lo	This rifle type projector is l .22 cal. Power loads. Multi	breech loaded and is powered by a .22 cal. Blank cartridge. ple ranges are achieved through the use of 4 different Powe ogth is 25 inches, the overall length is 41 inches and the we 5-75 yards. Size X-Long Range	er Loads. This model has a .50 cal.,
Projecto Arrow Ada	-	S to attach a Cap-Chur syringe to an arrow, for use with cross	s or longbows. Mfg. Cap-Chur
Blank Ada		with the outra long range prejector. It holds the power les	ad inside the breach

-	The Blank Adaptor is used with the extra-long range projector.	It holds the power-load inside the breech.
	KVS	Mfg.
	039-449	Cap-Chur

Carrying Case - Pistol

Zinnered closure. Convenient nockets for storing your syringes and equipment. Vinyl-covered nadded



KVSSizeMfg.039-463PistolCap-Chur	vinyi-covered, padded, carrying case. Zippered	i closure. Convenient pockets for	storing your synnges and equipment
039-463 Pistol Cap-Chur	KVS	Size	Mfg.
	039-463	Pistol	Cap-Chur



Guns and Accessories Projectors - Accessories **Carrying Case - Rifle** Vinyl-covered, padded, carrying case. Zippered closure. Convenient pockets for storing your syringes and equipment. **KVS** Size Mfg. 039-461 Rifle Cap-Chur **Plunger Lube** Premium silicone grease used to allow the rubber plunger to move forward freely and seal tightly along the barrel wall. SILGLYDE I **KVS** Size Mfg. 039-321 1.5oz Cap-Chur Positioner POSITIONER: is a clear fibreglass rod used to position the rubber plunger inside the syringe barrel. **KVS** Mfg. 039-445 Cap-Chur **Rubber Plunger** RUBBER PLUNGER: Comparable to the plunger on any syringe. Seals against the inside barrel wall and forces medication out through the needle. **KVS** Mfg. 039-443 Cap-Chur **Stock Extension** Shoulder stock for short range pistol. **KVS** Mfg. 039-309 Cap-Chur



Noseplugs

Noseplug Barb

NC NEEDLES: used primarily on thin skinned animals that are being "knocked out". Identified by the barb on the needle shank. Needle acts like a fish hook and will not come out in spite of animals activity.

777	KVS	Size	Part	Mfg.
	039-383	3/4"	NC1	Cap-Chur
,)	039-385	1 1/8"	NC2	Cap-Chur
	039-387	1 1/2"	NC3	Cap-Chur
	039-388	1 3/4"	NC4	Cap-Chur
	039-390	2"	NC5	Cap-Chur
	039-391	2.5"	NC6	Cap-Chur

Noseplug Heavy

Noseplugs (needles) for thick skinned animals. (500lbs & over). Diameter of the needle measures 0.125".

	KVS	Size	Part	Mfg.
	039-401	1 1/4"	NCL1	Cap-Chur
	039-403	1 1/2"	NCL2	Cap-Chur
7	039-405	1 3/4"	NCL3	Cap-Chur

Noseplug Large Collar

NS NEEDLE: Used for subcutaneous injections, identified by the short needle shank.

💧 К	(VS S	Size	Part	Mfg.
0	39-371	5/8"	NC1	Cap-Chur
0	39-373	3/4"	NS2	Cap-Chur

Noseplug Plain

PLAIN NEEDLES: Smooth shanked	d (no barb or collar) needle.	
KVS	Size	Mfg.
039-389	3/4"	Cap-Chur
039-393	1 1/8"	Cap-Chur

Noseplug Small Collar

NM NEEDLES: Used for IM injections, identified by the longer needle shank. The collar near the end of the NS & NM needles helps ensure the needle does not "pop out" of the animal before the injection is complete.

	KVS	Size	Part	Mfg.
	039-375	1 1/8"	NM1	Cap-Chur
	039-377	1 1/4"	NM2	Cap-Chur
	039-379	1 1/2"	NM3	Cap-Chur
	039-381	1 3/4"	NM4	Cap-Chur
_				•

Sideport

Sideport Noseplug - Drop Out

Sideport needle with drop out design.

KVS	Size	Mfg.
039-430	3/4"	Cap-Chur





Needles & Syringes

Sideport

Sideport Small Collar NEW: The side

NEW: The sideport delivery needle is designed primarily to eliminate the problem of "bounce outs" due to back pressure. Medication is dispensed over an area of about four times that of a regular needle.

	KVS	Size	Part	Mfg.
	039-425	1 1/8"	NM1	Cap-Chur
3	039-427	1 1/4"	NM2	Cap-Chur
	039-429	1 1/2"	NM3	Cap-Chur

Syringes

Marking Syringe Complete

For marking animals with colored ink using the Cap-Chur projectors as the delivery system. Marking syringes have a flat face noseplug, which does not pierce the skin of an animal. Upon impact with the animal, ink is squirted out of the exhaust ports in the syringe barrel walls, colouring the animals hide.

KVS	Size	Mfg.
039-475	15ml	Cap-Chur

Practice Syringe

syringe			
Practice syringe for the Cap-chur gun.			
KVS	Size	Mfg.	
039-447	1-3ml	Cap-Chur	

Syringe Barrel Only

The Cap-Chur syringe barrel is the holding tank or carrier of the injection while the syringe is in flight. The diameter of the syringe barrel must remain constant in order to fit through the projector. Therefore the length of the syringe barrel is proportional to how many ml's the barrel holds.

KVS	Size	Mfg.
039-347	1ml	Cap-Chur
039-349	2ml	Cap-Chur
039-351	3ml	Cap-Chur
039-353	4ml	Cap-Chur
039-355	5ml	Cap-Chur
039-357	7ml	Cap-Chur
039-359	10ml	Cap-Chur
039-361	15ml	Cap-Chur
039-363	20ml	Cap-Chur

Tail Piece

TAIL PIECE: is the "feather" or guide for the Cap-Chur syringe. Helps keep the syringe "true" in flight.

KVS	-	Mfg.	
039-441		Cap-Chur	



Powerloads & Charges

Charges

Charges

6

Charges are very small explosive devices which "fire" on impact. When the syringe hits the target, the rubber plunger is pushed forward by the expanding gas caused by the Cap Chur charge explosion. The drug is then forced out through the needle and into the target.

5 KVS	Size	Pack	Mfg.	TDG Class
039-451	1-3ml	50	Cap-Chur	LQ
039-453	4-10ml	50	Cap-Chur	LQ
039-455	15ml	50	Cap-Chur	LQ

Gas Packs

CO2 Gas Packs

Used in the long range + short range projectors to propel the Cap Chur syringes.

Brown

KVS	Pack	Mfg.
156-001	5	KVS

Powerloads

Powerloads - Low

POWER LOADS: .22 calibre blank loads used to propel Cap Chur syringes out of the Cap Chur Extra Long Range Projector. LOW: used at ranges from 15-40 yds.

ĸvs	Color	Size	Pack	Mfg.	TDG Class
039-313	Green	Low	50	Cap-Chur	LQ

Powerloads - Medium

POWER LOADS: .22 calibre blank loads used to propel Cap Chur syringes out of the Cap Chur Extra Long Range Projector. MEDIUM: used at ranges from 30-60 yds.

KVS	Color	Size	Pack	Mfg.	TDG Class
039-315	Yellow	Medium	50	Cap-Chur	LQ

Powerloads - Very Low

039-311

0

 POWER LOADS: .22 calibre blank loads used to propel Cap Chur syringes out of the Cap Chur Extra Long Range Projector.

 VERY LOW: used for ranges from 5-25 yds.

 KVS
 Color
 Size
 Pack
 Mfg.
 TDG Class

50

Cap-Chur

LQ

Very Low