

Cap-Chun

Notes



Table of Contents

Categorized Listing

| | |
|------------------------------|----|
| Guns and Accessories | 41 |
| Gun | 41 |
| Pistol | 41 |
| Rifle - Long Range | 41 |
| Rifle - X-Long Range | 41 |
| Gun Accessories | 41 |
| Arrow Adaptor | 41 |
| Blank Adaptor.22 | 41 |
| Carrying Case - Pistol | 41 |
| Carrying Case - Rifle | 41 |
| Plunger Lube | 42 |
| Positioner | 42 |
| Rubber Plunger | 42 |
| Stock Extension | 42 |
| Ultra Dot Scope | 42 |
| Needles & Syringes | 42 |
| Noseplugs | 42 |
| Noseplug Barb | 42 |
| Noseplug Heavy | 43 |
| Noseplug Large Collar | 43 |
| Noseplug Plain | 43 |
| Noseplug Small Collar | 43 |
| Sideport | 43 |
| Sideport Noseplug - Drop Out | 43 |
| Sideport Small Collar | 43 |
| Syringes | 44 |
| Marking Syringe Complete | 44 |
| Practice Syringe | 44 |
| Syringe Barrel Only | 44 |
| Tail Piece | 44 |
| Powerloads & Charges | 44 |
| Charges | 44 |
| Charges | 44 |
| Gas Packs | 45 |
| CO2 Gas Packs | 45 |
| Powerloads | 45 |
| Powerloads - High | 45 |
| Powerloads - Low | 45 |
| Powerloads - Medium | 45 |
| Powerloads - Very Low | 45 |

Table of Contents

Alphabetized Listing

| | | |
|------------------------------|--|-------|
| A | | |
| Arrow Adaptor | | 41 |
| B | | |
| Blank Adaptor.22 | | 41 |
| C | | |
| Carrying Case - Pistol | | 41 |
| Carrying Case - Rifle | | 41 |
| Charges | | 44 |
| CO2 Gas Packs | | 45 |
| G | | |
| Gas Packs | | 45 |
| Gun | | 41 |
| Gun Accessories | | 41-42 |
| Guns and Accessories | | 41-42 |
| M | | |
| Marking Syringe Complete | | 44 |
| N | | |
| Needles & Syringes | | 42-44 |
| Noseplug Barb | | 42 |
| Noseplug Heavy | | 43 |
| Noseplug Large Collar | | 43 |
| Noseplug Plain | | 43 |
| Noseplug Small Collar | | 43 |
| Noseplugs | | 42-43 |
| P | | |
| Pistol | | 41 |
| Plunger Lube | | 42 |
| Positioner | | 42 |
| Powerloads | | 45 |
| Powerloads - High | | 45 |
| Powerloads - Low | | 45 |
| Powerloads - Medium | | 45 |
| Powerloads - Very Low | | 45 |
| Powerloads & Charges | | 44-45 |
| Practice Syringe | | 44 |
| R | | |
| Rifle - Long Range | | 41 |
| Rifle - X-Long Range | | 41 |
| Rubber Plunger | | 42 |
| S | | |
| Sideport | | 43 |
| Sideport Noseplug - Drop Out | | 43 |
| Sideport Small Collar | | 43 |
| Stock Extension | | 42 |
| Syringe Barrel Only | | 44 |
| Syringes | | 44 |
| T | | |
| Tail Piece | | 44 |
| U | | |
| Ultra Dot Scope | | 42 |



Guns and Accessories

Gun

Pistol



Short range projector: a pistol type projector which is powered by one CO2 gaspack and can fire approx. 20 syringes for a distance of 40 ft. each. Mid range projector: a pistol type projector that is breech loaded and is powered by 2 - 12 gram CO2 gas packs.

| KVS | Size | Vendor |
|---------|-------------|----------|
| 039-305 | Short Range | Cap-Chur |
| 039-307 | Mid Range | Cap-Chur |

Rifle - Long Range



Long range projector: a rifle type projector which is powered by two CO2 gaspacks. The long range projector can be expected to fire 12 x 1cc syringes a distance of 25-30 yards with improved accuracy over the short range pistol.

| KVS | Vendor |
|---------|----------|
| 039-304 | Cap-Chur |

Rifle - X-Long Range



Extra long range projector: a rifle type projector which is powered by your choice of .22 blank power loads. The extra long range projector can be used at ranges from 5-75 yds (see power loads for more information on distances).

| KVS | Vendor | TDG Class |
|---------|----------|-----------|
| 039-301 | Cap-Chur | LQ |

Gun Accessories

Arrow Adaptor



This adaptor enables you to attach a Cap-Chur syringe to an arrow, for use with cross or longbows.

| KVS | Vendor |
|---------|----------|
| 039-465 | Cap-Chur |

Blank Adaptor.22



The Blank Adaptor is used with the extra-long range projector. It holds the power-load inside the breech.

| KVS | Vendor |
|---------|----------|
| 039-449 | Cap-Chur |

Carrying Case - Pistol



Vinyl-covered, padded, carrying case. Zippered closure. Convenient pockets for storing your syringes and equipment.

| KVS | Size | Vendor |
|---------|--------|----------|
| 039-463 | Pistol | Cap-Chur |

Carrying Case - Rifle



Vinyl-covered, padded, carrying case. Zippered closure. Convenient pockets for storing your syringes and equipment.

| KVS | Size | Vendor |
|---------|-------|----------|
| 039-461 | Rifle | Cap-Chur |



Guns and Accessories

Gun Accessories

Plunger Lube



Premium silicone grease used to allow the rubber plunger to move forward freely and seal tightly along the barrel wall.

| KVS | Size | Vendor |
|---------|-------|----------|
| 039-321 | 1.5oz | Cap-Chur |

Positioner



POSITIONER: is a clear fibreglass rod used to position the rubber plunger inside the syringe barrel.

| KVS | Vendor |
|---------|----------|
| 039-445 | Cap-Chur |

Rubber Plunger



RUBBER PLUNGER: Comparable to the plunger on any syringe. Seals against the inside barrel wall and forces medication out through the needle.

| KVS | Vendor |
|---------|----------|
| 039-443 | Cap-Chur |

Stock Extension



Shoulder stock for short range pistol.

| KVS | Vendor |
|---------|----------|
| 039-309 | Cap-Chur |

Ultra Dot Scope



Laser sight for Cap-Chur guns.

| KVS | Vendor |
|---------|----------|
| 039-319 | Cap-Chur |

Needles & Syringes

Noseplugs

Noseplug Barb



NC NEEDLES: used primarily on thin skinned animals that are being "knocked out". Identified by the barb on the needle shank. Needle acts like a fish hook and will not come out in spite of animals activity.

| KVS | Size | Part | Vendor |
|---------|-----------|------|----------|
| 039-383 | 3/4inch | NC1 | Cap-Chur |
| 039-385 | 1 1/8inch | NC2 | Cap-Chur |
| 039-387 | 1 1/2inch | NC3 | Cap-Chur |
| 039-388 | 1 3/4inch | NC4 | Cap-Chur |
| 039-390 | 2inch | NC5 | Cap-Chur |
| 039-391 | 2.5inch | NC6 | Cap-Chur |



Needles & Syringes

Noseplugs

Noseplug Heavy

Noseplugs (needles) for thick skinned animals. (500lbs & over). Diameter of the needle measures 0.125".



| KVS | Size | Part | Vendor |
|---------|-----------|------|----------|
| 039-401 | 1 1/4inch | NCL1 | Cap-Chur |
| 039-403 | 1 1/2inch | NCL2 | Cap-Chur |
| 039-405 | 1 3/4inch | NCL3 | Cap-Chur |

Noseplug Large Collar

NS NEEDLE: Used for subcutaneous injections, identified by the short needle shank.



| KVS | Size | Part | Vendor |
|---------|---------|------|----------|
| 039-371 | 5/8inch | NC1 | Cap-Chur |
| 039-373 | 3/4inch | NS2 | Cap-Chur |

Noseplug Plain

PLAIN NEEDLES: Smooth shanked (no barb or collar) needle.



| KVS | Size | Vendor |
|---------|-----------|----------|
| 039-389 | 3/4inch | Cap-Chur |
| 039-393 | 1 1/8inch | Cap-Chur |

Noseplug Small Collar

NM NEEDLES: Used for IM injections, identified by the longer needle shank. The collar near the end of the NS & NM needles helps ensure the needle does not "pop out" of the animal before the injection is complete.



| KVS | Size | Part | Vendor |
|---------|-----------|------|----------|
| 039-375 | 1 1/8inch | NM1 | Cap-Chur |
| 039-377 | 1 1/4inch | NM2 | Cap-Chur |
| 039-379 | 1 1/2inch | NM3 | Cap-Chur |
| 039-381 | 1 3/4inch | NM4 | Cap-Chur |

Sideport

Sideport Noseplug - Drop Out

Sideport needle with drop out design.



| KVS | Size | Vendor |
|---------|------|----------|
| 039-430 | 3/4" | Cap-Chur |

Sideport Small Collar

NEW: The sideport delivery needle is designed primarily to eliminate the problem of "bounce outs" due to back pressure. Medication is dispensed over an area of about four times that of a regular needle.



| KVS | Size | Part | Vendor |
|---------|-----------|------|----------|
| 039-425 | 1 1/8inch | NM1 | Cap-Chur |
| 039-427 | 1 1/4inch | NM2 | Cap-Chur |
| 039-429 | 1 1/2inch | NM3 | Cap-Chur |



Needles & Syringes

Syringes

Marking Syringe Complete



For marking animals with colored ink using the Cap-Chur projectors as the delivery system. Marking syringes have a flat face noseplug, which does not pierce the skin of an animal. Upon impact with the animal, ink is squirted out of the exhaust ports in the syringe barrel walls, colouring the animals hide.

| KVS | Size | Vendor |
|---------|------|----------|
| 039-475 | 15ml | Cap-Chur |

Practice Syringe



Practice syringe for the Cap-chur gun.

| KVS | Size | Vendor |
|---------|-------|----------|
| 039-447 | 1-3ml | Cap-Chur |

Syringe Barrel Only



The Cap-Chur syringe barrel is the holding tank or carrier of the injection while the syringe is in flight. The diameter of the syringe barrel must remain constant in order to fit through the projector. Therefore the length of the syringe barrel is proportional to how many ml's the barrel holds.

| KVS | Size | Vendor |
|---------|------|----------|
| 039-347 | 1ml | Cap-Chur |
| 039-349 | 2ml | Cap-Chur |
| 039-351 | 3ml | Cap-Chur |
| 039-353 | 4ml | Cap-Chur |
| 039-355 | 5ml | Cap-Chur |
| 039-357 | 7ml | Cap-Chur |
| 039-359 | 10ml | Cap-Chur |
| 039-361 | 15ml | Cap-Chur |
| 039-363 | 20ml | Cap-Chur |

Tail Piece



TAIL PIECE: is the "feather" or guide for the Cap-Chur syringe. Helps keep the syringe "true" in flight.

| KVS | Vendor |
|---------|----------|
| 039-441 | Cap-Chur |

Powerloads & Charges

Charges

Charges



Charges are very small explosive devices which "fire" on impact. When the syringe hits the target, the rubber plunger is pushed forward by the expanding gas caused by the Cap Chur charge explosion. The drug is then forced out through the needle and into the target.

| KVS | Size | Pack | Vendor | TDG Class |
|---------|--------|------|----------|-----------|
| 039-451 | 1-3ml | 50 | Cap-Chur | LQ |
| 039-453 | 4-10ml | 50 | Cap-Chur | LQ |
| 039-455 | 15ml | 50 | Cap-Chur | LQ |



Powerloads & Charges Gas Packs

CO2 Gas Packs



Used in the long range + short range projectors to propel the Cap Chur syringes.

| KVS | Pack | Vendor |
|---------|------|--------|
| 156-001 | 5 | KVS |

Powerloads

Powerloads - High



POWER LOADS: .22 calibre blank loads used to propel Cap Chur syringes out of the Cap Chur Extra Long Range Projector. HIGH : used at ranges from 40-75 yds.

| KVS | Color | Size | Pack | Vendor | TDG Class |
|---------|-------|------|------|----------|-----------|
| 039-317 | Red | High | 50 | Cap-Chur | LQ |

Powerloads - Low



POWER LOADS: .22 calibre blank loads used to propel Cap Chur syringes out of the Cap Chur Extra Long Range Projector. LOW: used at ranges from 15-40 yds.

| KVS | Color | Size | Pack | Vendor | TDG Class |
|---------|-------|------|------|----------|-----------|
| 039-313 | Green | Low | 50 | Cap-Chur | LQ |

Powerloads - Medium



POWER LOADS: .22 calibre blank loads used to propel Cap Chur syringes out of the Cap Chur Extra Long Range Projector. MEDIUM: used at ranges from 30-60 yds.

| KVS | Color | Size | Pack | Vendor | TDG Class |
|---------|--------|--------|------|----------|-----------|
| 039-315 | Yellow | Medium | 50 | Cap-Chur | LQ |

Powerloads - Very Low



POWER LOADS: .22 calibre blank loads used to propel Cap Chur syringes out of the Cap Chur Extra Long Range Projector. VERY LOW: used for ranges from 5-25 yds.

| KVS | Color | Size | Pack | Vendor | TDG Class |
|---------|-------|----------|------|----------|-----------|
| 039-311 | Brown | Very Low | 50 | Cap-Chur | LQ |

CAPCHUR

